ELISA bLF kit

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**ELISA kit for quantification of bovine lactoferrin in
milk and milk's powders, colostrum and
colostrum's powders and in bovine sera**

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Kit ELISA de dosage de la lactoferrine bovine

Catalogue reference E_bLF_01

ELISA kit for quantification of bovine lactoferrin

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**ELISA kit for quantification of bovine lactoferrin in
milk and milk's powders, colostrum and colostrum's
powders and in bovine sera**

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**This document must be fully read before running the kit for
the first time.**

Quality Management System **ISO9001:2008** certified by

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TABLE OF CONTENTS

1. INTRODUCTION	17
2. FIELD OF APPLICATION	17
3. ASSAY PRINCIPLE	17
4. KIT CONTENT.....	19
5. STORAGE CONDITIONS	19
6. REAGENTS AND EQUIPMENTS REQUIRED	20
7. PREPARATION OF THE WASH BUFFER	20
8. SAMPLES PREPARATION	21
9. REAGENTS PREPARATION.....	23
9.1. STANDARDS PREPARATION.....	23
9.2. CONJUGATE SOLUTION PREPARATION (IgG-HRP).....	23
10. INSTRUCTIONS.....	24
11. RESULTS.....	27
12. TECHNICAL ADVICES	28
13. SAFETY RECOMMENDATIONS	28
14. ASSAY CHARACTERISTICS.....	29
ANNEXE / ANNEX.....	30

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1. INTRODUCTION

Bovine lactoferrin (**bLF**) is a 80 kDa glycoprotein able to bind iron (Fe^{3+}). It is found in most of the biological fluids, mainly in milk, and plays key roles in the defense mechanism of the organism. A lot of biological activities have been reported for **bLF**: antimicrobial, immuno-modulator, anti-tumoral, absorption of iron or bone growth. Antioxidant, antiviral, antifungal and anti-inflammatory activities have also been reported for **bLF**.

In November 2012, **bLF** purified from milk and whey was approved as a « Novel Food Ingredient » by the European Food Safety Authority (EFSA). Accordingly **bLF** can now be used as an additive in many food products such as infant formula and dairy products.

2. FIELD OF APPLICATION

This kit is intended for the quantification of bovine lactoferrin (bLF) in:

- bovine milk and milk powders
- bovine colostrum and colostrum powder
- bovine serum


This kit is intended for research and manufacturing use only, not for diagnostic uses.

3. ASSAY PRINCIPLE

This immuno-enzymatic assay is a sandwich ELISA.

It was developed from a rabbit serum raised against bovine lactoferrin (bLF). The immunoglobulins (IgG) were purified by ion exchange chromatography. A fraction of these IgG is used as purified capture antibody: the ELISA plate is coated with the purified fraction of IgG, stabilized and delivered ready-to-use.

Another fraction was labeled to horseradish peroxidase (IgG-HRP) and is used as a conjugate antibody.

When performing the assay, wells of the ELISA plate are filled with standards and samples. The anti-bLF antibody solid phase capture the bovine lactoferrin contained in the standards and samples. The ELISA plate is washed to remove any unbound material.

The solution of IgG-HRP is distributed in ELISA wells and react with the antibody-bLF complex of the solid phase. The reaction results in the formation of a sandwich complex. After washing, the complex is revealed using tetramethyl benzidine chromogen substrate (TMB).

The Optical Density (OD) generated by the hydrolysis of the substrate is measured using a microtiter plate reader. The OD is directly proportional to the concentration of bovine lactoferrin .

4. KIT CONTENT

The kit is delivered in a box containing:

- 6 strips of 16wells coated with purified anti-bLF immunoglobulins in a foil, zip-lock bag with dessicant.
NB : a preliminary wash of the ELISA plate is necessary before performing the assay (cf § 10.)
- **R1** : Wash buffer 10-fold concentrated – 1 x 100mL vial
- **R2** : Dilution buffer for standards, samples and the conjugate, ready to use – 4 x 30mL vial
- **R3** : Standard solution at 20µg/ml of bovine lactoferrin – 1 x 250µL microtube
- **R4** : conjugate solution, purified IgG labeled to HRP (Ig-HRP), 40-fold concentrated – 1 x 500µL microtube
- **SB** : tetramethyl benzidine substrate (TMB), ready-to-use – 1 x 15mL vial
- **ST** : stop solution H₂SO₄ 0,45M, ready to use – 1 x 10mL vial
- 1 adhesive

5. STORAGE CONDITIONS

Store all reagents between +2°C to +8°C.

Do not use the reagents after their expiry date.

Prevent reagents R4 and SB from light exposure during storage and assay incubation steps.

Before use, bring all reagents to room temperature.

Gently homogenise each vial manually by turning upside down several times before use.

Do not mix reagents from different lots.

6. REAGENTS AND EQUIPMENTS REQUIRED

- Distilled or deionised water (PANREAC 212236.1214, or equivalent)
- Disposable glass tubes
- Disposable 1.2ml Ø8.8mm x 45mm microtitration tubes with rack
- Container to prepare wash buffer
- Reservoirs for the delivery of the reagents (for multichannel pipette)
- Adhesive films
- Precision micropipettes , 10 – 5000 µl with adapted disposable tips
- 8-channel micropipette with variable volume (50-200 µL) with adapted disposable tips
- Microtiter plate washer (manual or automatic)
- Orbital shaker (type IKA-Vibrax® or equivalent) with tray for microtiter plate
- Microtiter plate reader with 450nm filter (Biotek ELx808 or equivalent)
- Data treatment software including the 4 parameter logistics regression model $Y = (A-D)/(1+(X/C)^B) + D$

7. PREPARATION OF THE WASH BUFFER

Prepare the wash buffer as follows:

Into a clean vial, mix the 100mL of 10X wash buffer **R1** with 900mL of distilled or deionised water.

1. Homogenise slowly
2. This solution can be stored at +2°C/+8°C for two weeks

8. SAMPLES PREPARATION

- If samples will not be assayed immediately, stored refrigerated for up to a few days or frozen for long-term storage (avoid freeze-thaw cycles)
- Bring to room temperature and vortex the samples before preparation.
- **Samples may be assayed with proper dilution. Dilute the samples, based on the expected concentration of bovine lactoferrin, to fall within the range of the standard curve.**

sample	Concentration bLF ($\mu\text{g/ml}$)	Recommended dilution
raw milk	150 – 200	1/4000
pasteurized milk crude / semi-skimmed	50 – 100	1/2000
sterilised milk	< 0,03	1/10
colostrum	700 – 900	1/20000
normal bovine serum	0,5 – 0,7	1/20
milk powder whole and semi-skimmed reconstituted at 100mg/mL	0,5 – 2	1/25
milk powder skimmed reconstituted at 100 mg/mL	10-30	1/400
colostrum powder reconstituted at 100 mg/mL	200 – 300	1/4000

The recommended dilution for samples should be used as a guideline. If unsure of sample level, a serial dilution with one or two representative samples before running the entire plate is highly recommended.

In order to optimize the use of the dilution buffer R2 provided in the kit, we advise to dilute the samples by serial dilution under appropriate volumes.

For example, a dilution 1/4000 can be realized as follow:

dil1	1/10	50 μl (sample)	+ 450 μl buffer R2
dil2	1/10 (1/100 final)	50 μl (dil1)	+ 450 μl buffer R2
dil3	1/10 (1/1000 final)	50 μl (dil2)	+ 450 μl buffer R2
dil4	1/4 (1/4000 final)	125 μl (dil3)	+ 375 μl buffer R2

Reconstitute milk and colostrum powders at 100 mg/mL in deionised water. Homogenise with a vortex and mix slowly on a rotative shaker for 20min at room temperature. Centrifuge 5min (>10000g), collect the supernatant and then dilute at the recommended dilution.

It is recommended that each laboratory uses internal quality controls samples in routine analysis to insure that all reagents and procedures are correct.

9. REAGENTS PREPARATION

9.1. Standards preparation

A working solution at 2000 ng/mL of bovine lactoferrin is prepared extemporaneously by diluting the stock solution **R3** (solution at 20µg/ml of bovine lactoferrin) to 1/10 in the **R2** dilution buffer as follows:

50µL (stock solution **R3**) + 450µL **R2** buffer

Homogenise the solution.

Discard after use.

The standard solutions are prepared in micro-titration tubes as follows:

STD 1 => 200 ng/ml	: 100 µL (work. Sol.) +	900 µL R2 buffer
STD 2 => 100 ng/ml	: 500 µL (STD 1) +	500 µL R2 buffer
STD 3 => 50 ng/ml	: 500 µL (STD 2) +	500 µL R2 buffer
STD 4 => 25 ng/ml	: 500 µL (STD 3) +	500 µL R2 buffer
STD 5 => 12,5 ng/ml	: 500 µL (STD 4) +	500 µL R2 buffer
STD 6 => 6,25 ng/ml	: 500 µL (STD 5) +	500 µL R2 buffer
STD 7 => 3,125 ng/ml	: 500 µL (STD 6) +	500 µL R2 buffer
Blank (BLK)	:	500 µL R2 buffer

Homogenise the solution with an orbital shaker between each dilution.

Discard after use.

The volumes are indicated for standards assays in duplicate .

9.2. Conjugate solution preparation (Ig-HRP)

Prepare a 1/40 diluted solution of the 40-fold concentrated Conjugate Ig-HRP **R4** in the dilution buffer **R2**.

Ig_HRP 1/40: 300µL (**R4** conjugate 40X) + 11,7mL **R2** buffer

Homogenise slowly by manual inversions of the vial.

Discard after use.

10. INSTRUCTIONS

Bring all reagents to room temperature, including the coated ELISA plate in its aluminum pouch kept closed (do not open until temperature equilibration). All standards, controls, and samples should be tested at least in duplicate.

1. Prepare the wash buffer, samples, standards and controls. Remove the appropriate number of strips from the aluminum pouch and return to +2°C/+8°C. Take care not to touch the underside of the plate. Place all the microtitration tubes containing the standards, blanks, samples and controls on a rack following the same location of the plate layout.

2. Wash the ELISA strips with an automated or manual plate washer as following:

Distribute 250µl to 300µl of 1X wash buffer (R1) in each well. Drain the wells by dumping the plate or aspiration system. Repeat for a total of 4 washes. After the final wash, tap over absorbent paper to remove most of the residual liquid.

NB : Do not allow wells to dry between each step of the assay.

3. Incubation of standards, blanks, samples and controls.

Pipette **100µL/well** (use an 8-channel pipette) of the standards, blanks, samples and controls.

Cover with an adhesive and shake gently for a few seconds on a plate shaker (<500 rpm).

Incubate 1 hour at 25°C without agitation.

Wash 5 times as described in section 2.

4. Incubation of the solution of conjugate (IgG-HRP)

Prevent the solution and filled plate from light exposure during this step.

Pipette **100µL/well** of the conjugate solution (use an 8-channel pipette).

Cover with an adhesive and shake gently for few seconds on a plate shaker (<500 rpm).

Incubate 1 hour at 25°C without agitation in the dark.

Wash 5 times as described in section 2.

5. Incubation of the substrate (TMB)

Discard if the substrate has a distinct blue color prior to the assay.

Pipette **100µL/well** of the ready-to-use substrate solution (use an 8-channel pipette).

Cover with an adhesive and shake gently for a few seconds on a plate shaker (<500 rpm).

Incubate 30minutes at 25°C without agitation in the dark.

A blue coloration will develop gradually during this step.

6. Stop solution

Pour the ready-to-use stop solution **ST** H_2SO_4 0,45M into a clean reservoir.

Pipette **50 μl /well** of the stop solution (use an 8-channel pipette). The blue coloration will turn to yellow.

Cover with an adhesive and shake gently for a few seconds on a plate shaker (<500 rpm).

7. Remove the adhesive and read the absorbance at 450nm with a microplate reader.

Simplified assay protocol

Prepare the wash buffer, standards, samples and controls

Wash the strips

Pipette **100µl/well** of standards, blanks,
samples and controls in duplicate
Incubate **1 hour** at 25°C without agitation

Wash the strips

Pipette **100µl/well** of the conjugate solution (IgG-HRP)
Incubate **1 hour** at 25°C without agitation in the dark

Wash the strips

Pipette **100µl/well** of the substrate
Incubate from **15 to 30 minutes** at 25°C in the dark

Fill **50µl** of the stop solution


Read the absorbance at 450nm

11. RESULTS

The standard curve is constructed by reporting the OD at 450nm of each standard on the Y-axis versus their respective concentrations of HCP in ng/mL on the X-axis. The data reduction is performed through a computer software using curve fitting models such as 4 parameter logistic model ($Y = (A-D)/(1+(X/C)^B)+D$).

Exemple of standard curve:

Absorbances of samples are interpolated from the calibration curve.

It is recommended that each laboratory uses internal quality controls samples in each routine analysis to insure that all reagents and procedures are correct.

12. TECHNICAL ADVICES

Any improper handling of samples or modification of the test procedure may influence the results. The indicated pipetting volumes, incubation times, temperatures and pretreatment steps have to be performed strictly according to the instructions.

Avoid contamination of reagents, pipettes and wells/tubes. Use new disposable plastic pipette tips for each component and sample. Do not interchange caps. Always cap not used vials. Do not reuse wells/tubes or reagents.

Once the test has been started, all steps should be completed without interruption. Make sure that required reagents, materials and devices are prepared ready at the appropriate time. Allow all reagents and samples to reach room temperature and gently homogenize each vial of liquid reagent and sample before use. Mix reagents without foaming.

Incubation time affects results. All wells should be handled in the same order and time sequences. It is recommended to use an 8-channel Micropipettor for pipetting of solutions in all wells. **Maintain a repetitive timing sequence from well to well for all assay steps to insure that all incubation times are the same for each well.**

Microplate washing is important. Improperly washed wells will give erroneous results. It is recommended to use a multichannel pipette or an automatic microplate washing system. Do not allow the wells to dry between incubations. Do not scratch coated wells during rinsing and aspiration. Rinse and fill all reagents with care. While rinsing, check that all wells are filled precisely with Wash Buffer, and that there are no residues in the wells.

Use qualified and approved equipment and devices (pipettes, plate reader,...).

Humidity can affect the ELISA coated wells. Do not open the aluminum pouch until the ELISA plate reaches room temperature.

Change pipette tips between each sample, calibration standard and reagent to avoid cross-contamination.

13. SAFETY RECOMMENDATIONS

Follow good laboratory practice and safety guidelines.

Wear lab coats, disposable gloves and protective glasses where/when necessary.

Stop solution is 0,45M H₂SO₄. Avoid contact with eyes, skin, and clothing. At the concentrations used in this kit, none of the other reagents are believed to be harmful.

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Linearity	non linear regression 4 parameters logistic model $r^2 > 0,99$
Limit of detection	LOD = 1,8 ng/ml interpolation of mean absorbance of blank + 3 sd
Limit of quantification	LOQ = 3,0 ng/ml interpolation of mean absorbance of blank + 10 sd
Repeatability (intra-assay)	CV% = 1 - 5% at concentration levels 25, 50 and 100 ng/mL for milk, colostrum powder and serum within run n=6, 3 separate run-to-run assays
Intermediate precision (inter-assay)	CV% = 5 - 14% at concentration levels 25, 50 and 100 ng/mL for milk, colostrum powder and serum within run n=6, 3 separate run-to-run assays
Exactitude	% of recovery = 95 - 140% for milk (115% ± 20%) Spike at 3 concentrations level in milk
Range of quantification	Range of quantification : 5 – 100 ng/ml
Specificity	Cross reaction : weak with goat or sheep milk and sow colostrum medium with buffalo milk no cross reaction with camel milk and mare colostrum No cross reaction with bovine whey protein (betalactoglobulin, alphasalactalbumin, immunoglobulin and albumine)

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	1	2	3	4	5	6	7	8	9	10	11	12
A												
B												
C												
D												
E												
F												
G												
H												

IDBiotech

*ImmunoDiffusion Biotechnologies SARL
Avenue Marie Curie
63500 Issoire – FRANCE*

Tel: +33 (0)4 73 54 95 01

Fax: +33 (0)4 73 54 43 94

Web : www.idbiotech.com

e-mail : contact@idbiotech.com